


Anti-Pdk3 Polyclonal Antibody For Wb, Elisa - Q15120

Item Number: CM0018537



Introduction

Rabbit polyclonal antibody against human PDK3, validated for Western blot and ELISA. Targets mitochondrial pyruvate dehydrogenase kinase isozyme 3, a regulator of glucose metabolism and aerobic respiration.

[Learn More](#)

Parameter	Value
Product Name	PDK3 Rabbit pAb
Remarks/Alias	CMTX6; GS1-358P8.4; PDK3
Species	Human
GeneID (Human)	5165
GeneID	5165
Immunogen	Recombinant protein, Recombinant fusion protein containing a sequence corresponding to amino acids 137-406 of human PDK3 (NP_005382.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q15120
Protein Weight	47kDa
Shipping	Ice bag