


Anti-Pdlim1/Clp36 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - O00151

Item Number: CM0003161



Introduction

Rabbit monoclonal PDLIM1/CLP36 antibody, validated for WB, IF/ICC, and ELISA. Cross-reacts with human and mouse. Ideal for cytoskeletal adapter studies, stress fiber dynamics, and cell migration research.

[Learn More](#)

Parameter	Value
Product Name	PDLIM1/CLP36 Rabbit mAb
Remarks/Alias	CLIM1; CLP36; CLP-36; hCLIM1; HEL-S-112; PDLIM1/CLP36
Species	Human
GeneID (Human)	9124
GeneID	9124
Immunogen	Synthetic peptide; A synthetic peptide corresponding to a sequence within amino acids 200-300 of human PDLIM1/CLP36 (O00151)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	O00151
Protein Weight	36kDa
Shipping	Ice bag