


Anti-Pdlim4 Polyclonal Antibody For Wb, If/Icc, Elisa - P50479

Item Number: CM0022657



Introduction

PDLIM4 Rabbit pAb (P50479) is a polyclonal antibody validated for WB, IF/ICC, and ELISA in human, mouse, and rat. Targets PDZ and LIM domain protein 4, which suppresses SRC activation and regulates actin cytoskeleton reorganization.

[Learn More](#)

Parameter	Value
Product Name	PDLIM4 Rabbit pAb
Remarks/Alias	RIL; PDLIM4
Species	Human
GeneID (human)	8572
GeneID	8572
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-210 of human PDLIM4 (NP_003678.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50479
Protein Weight	35kDa
Shipping	Ice bag