

## Anti-Pdlim5 Rabbit Polyclonal Antibody For Wb And Elisa - Q96Hc4

Item Number: CM0030173



### Introduction

Rabbit polyclonal antibody targeting human PDLIM5, suitable for WB and ELISA. Cross-reacts with mouse, rat. Ideal for studying heart development, cardiomyopathy, and neuropsychiatric disorders.

[Learn More](#)

| Parameter        | Specification                                                                                                       |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | PDLIM5 Rabbit pAb                                                                                                   |
| Remarks/Alias    | L9; ENH; LIM; ENH1; PDLIM5                                                                                          |
| Species          | Human                                                                                                               |
| GeneID (Human)   | 10611                                                                                                               |
| GeneID           | 10611                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 317-596 of human PDLIM5 (NP_006448.4) |
| Source           | Rabbit                                                                                                              |
| Category         | Polyclonal Antibodies                                                                                               |
| Application      | WB, ELISA                                                                                                           |
| Cross Reactivity | Human, Mouse, Rat                                                                                                   |
| SWISS            | Q96HC4                                                                                                              |
| Protein Weight   | 64kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |