

Anti-Pdlim5 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Hc4

Item Number: CM0019234



Introduction

Rabbit polyclonal PDLIM5 antibody validated for WB and ELISA. Reacts with human. PDLIM5 scaffolds PKC at cardiac Z-disks and shapes dendritic spines; implicated in bipolar disorder, schizophrenia. Ideal for heart and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	PDLIM5 Rabbit pAb
Remarks/Alias	L9; ENH; LIM; ENH1; PDLIM5
Species	Human
Gene ID	10611
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 300-400 of human PDLIM5 (NP_006448.5).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96HC4
Protein Weight	64kDa
Shipping	Ice bag