

Anti-Pdp1 Polyclonal Antibody For Wb - Q9P0J1

Item Number: CM0015703

Introduction

Rabbit polyclonal antibody against human PDP1, a mitochondrial pyruvate dehydrogenase phosphatase. Suitable for WB, IF/ICC, ELISA; cross-reacts with human, mouse, rat. Predicted MW 61 kDa. Ideal for metabolism and mitochondrial research.

[Learn More](#)

Parameter	Value
Product Name	PDP1 Rabbit pAb
Remarks/Alias	PDH; PDP; PDPC; PPM2A; PPM2C; PDP1
Species	Human
GeneID (human)	54704
GeneID	54704
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 72-260 of human PDP1 (NP_060914.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9P0J1
Protein Weight	61kDa
Shipping	Ice bag