

Anti-Pdp2 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9P2J9

Item Number: CM0024591

Introduction

High-quality rabbit polyclonal antibody against human PDP2, validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for studying mitochondrial pyruvate dehydrogenase regulation and immunometabolism.

[Learn More](#)

Parameter	Value
Product Name	PDP2 Rabbit pAb
Remarks/Alias	PPM2B, PDPC 2, PPM2C2, PDP2
Species	Human
GeneID (Human)	57546
GeneID	57546
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 67-250 of human PDP2 (NP_065837.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9P2J9
Protein Weight	60kDa
Shipping	Ice bag