

Anti-Pdx1 Polyclonal Antibody For Wb And Elisa - P52945

Item Number: CM0014294



Introduction

Rabbit anti-PDX1 polyclonal antibody validated for Western blot and ELISA, recognizing human and mouse PDX1 (31 kDa). Ideal for diabetes and pancreatic development research.

[Learn More](#)

Parameter	Value
Product Name	PDX1 Rabbit pAb
Remarks/Alias	GSF, IPF1, IUF1, IDX-1, MODY4, PDX-1, STF-1, PAGEN1, PDX1
Species	Human
Gene ID	3651
Immunogen	Synthetic peptide corresponding to amino acids 1-100 of human PDX1 (NP_000200.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P52945
Protein Weight	31kDa
Shipping	Ice bag