


Anti-Pdxk Monoclonal Antibody For Wb, Elisa - 000764

Item Number: CM0007272



Introduction

Recombinant anti-PDXK rabbit monoclonal antibody validated for WB and ELISA. Detects human PDXK (35 kDa, SWISS O00764). Ideal for vitamin B6 metabolism and kinase research.

[Learn More](#)

Parameter	Value
Product Name	PDXK Rabbit mAb
Remarks/Alias	PKH, PNK, HMSN6C, PRED79, C21orf97, HEL-S-1a, C21orf124, PDXK
Species	Human
GeneID (Human)	8566
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-312 of human PDXK (NP_003672.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O00764
Protein Weight	35kDa
Shipping	Ice bag