

Anti-Pdxk Polyclonal Antibody For Wb, If/Icc, Elisa - 000764

Item Number: CM0029155



Introduction

High-quality rabbit polyclonal anti-PDXK antibody validated for WB, IF/ICC, ELISA. Detects human, mouse, rat PDXK. Ideal for vitamin B6 metabolism, kinase activity, and neuropathy research.

[Learn More](#)

Parameter	Value
Product Name	PDXK Rabbit pAb
Remarks/Alias	PKH; PNK; HMSN6C; PRED79; C21orf97; HEL-S-1a; C21orf124; PDXK
Species	Human
Gene ID	8566 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-312 of human PDXK (NP_003672.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00764
Protein Weight	35kDa
Shipping	Ice bag