

Anti-Pdxd Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q96Gd0

Item Number: CM0020146



Introduction

Rabbit polyclonal antibody against PDXP (Chronophin) for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat PDXP (32 kDa). Useful for studying vitamin B6 metabolism, cofilin signaling, and neuronal development.

[Learn More](#)

Parameter	Value
Product Name	PDXP Rabbit pAb
Remarks/Alias	CIN; PLP; dJ37E16.5; PDXP
Species	Human
Gene ID (Human)	57026
Gene ID	57026
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 200 to the C-terminus of human PDXP (NP_064711.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96GD0
Protein Weight	32kDa
Shipping	Ice bag