

Anti-Pdzk1 Rabbit Polyclonal Antibody For Wb, Elisa - Q5T2W1

Item Number: CM0025579



Introduction

Rabbit polyclonal antibody against human PDZK1 (NHERF3/CAP70) scaffold protein, validated for Western blot and ELISA, cross-reactive with mouse. Ideal for studying epithelial ion transport and membrane protein regulation. Offered by CamelBio.

[Learn More](#)

Parameter	Value
Product Name	PDZK1 Rabbit pAb
Remarks/Alias	CAP70; CLAMP; PDZD1; NHERF3; NHERF-3; PDZK1
Species	Human
GeneID (human)	5174
GeneID	5174
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human PDZK1 (NP_001188254.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q5T2W1
Protein Weight	57kDa
Shipping	Ice bag