


## Anti-Pebp1 Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P30086

Item Number: CM0016704



### Introduction

Anti-PEBP1 rabbit polyclonal antibody for WB, IF/ICC, IP, ELISA. Targets human PEBP1 (RKIP) with cross-reactivity to mouse and rat. Ideal for kinase signaling and CNS research.

[Learn More](#)

| Parameter        | Value                                                                                                                                 |
|------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | PEBP1 Rabbit pAb                                                                                                                      |
| Remarks/Alias    | PBP; HCNP; PEBP; RKIP; HCNPPp; PEBP-1; HEL-210; HEL-S-34; HEL-S-96; PEBP1                                                             |
| Species          | Human                                                                                                                                 |
| Gene ID (Human)  | 5037                                                                                                                                  |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-187 of human PEBP1 (NP_002558.1). |
| Source           | Rabbit                                                                                                                                |
| Category         | Polyclonal Antibodies                                                                                                                 |
| Application      | WB, IF/ICC, IP, ELISA                                                                                                                 |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                     |
| SWISS            | P30086                                                                                                                                |
| Protein Weight   | 21kDa                                                                                                                                 |
| Shipping         | Ice bag                                                                                                                               |