


Anti-Pebp1 Rabbit Monoclonal Antibody For Wb, Elisa - P30086

Item Number: CM0002649



Introduction

Rabbit monoclonal antibody against human PEBP1 (RKIP/HCNP). Validated for WB and ELISA, cross-reacts with mouse and rat. Ideal for MAPK signaling and neurobiology research.

[Learn More](#)

Parameter	Value
Product Name	PEBP1 Rabbit mAb
Remarks/Alias	PBP; HCNP; PEBP; RKIP; HCNPpp; PEBP-1; HEL-210; HEL-S-34; HEL-S-96; PEBP1
Species	Human
Gene ID (Human)	5037
Gene ID	5037
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PEBP1 (P30086).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30086
Protein Weight	21kDa
Shipping	Ice bag