

Anti-Pebp4 Polyclonal Antibody For Wb And Elisa - Q96S96

Item Number: CM0016739

Introduction

Rabbit polyclonal antibody against human PEBP4 (CORK1), validated for WB and ELISA, with mouse cross-reactivity. Ideal for PI3K-AKT signaling and secreted protein research.

[Learn More](#)

Parameter	Specification
Product Name	PEBP4 Rabbit pAb
Remarks/Alias	CORK1; CORK-1; PEBP-4; hPEBP4; PRO4408; GWTM1933; HEL-S-300; PEBP4
Species	Human
GeneID (Human)	157310
GeneID	157310
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q96S96
Protein Weight	26kDa
Shipping	Ice bag