

Anti-Pebp4 Rabbit Polyclonal Antibody - Q9D9G2

Item Number: CM0021977

Introduction

Rabbit polyclonal antibody against PEBP4 (Phosphatidylethanolamine-binding protein 4), validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for PI3K-AKT signaling pathway research.

[Learn More](#)

Parameter	Value
Product Name	PEBP4 Rabbit pAb
Remarks/Alias	PEBP-4; PEBP4
Species	Mouse
GeneID (human)	157310
GeneID	73523
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 90-242 of mouse PEBP4 (NP_082836.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9D9G2
Protein Weight	27kDa
Shipping	Ice bag