


Anti-Pedf/Serpinf1 Polyclonal Antibody For Wb, If/Icc, Elisa - P36955

Item Number: CM0019882



Introduction

Rabbit polyclonal antibody against human PEDF/SERPINF1 (P36955). Validated for WB, IF/ICC, ELISA. Cross-reacts with mouse. Ideal for angiogenesis and neurotrophic research.

[Learn More](#)

Parameter	Value
Product Name	PEDF/SERPINF1 Rabbit pAb
Remarks/Alias	OI6; OI12; PEDF; EPC-1; PIG35
Species	Human
Gene ID (Human)	5176
Gene ID	5176
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 141-240 of human PEDF/SERPINF1 (NP_002606.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P36955
Protein Weight	46kDa
Shipping	Ice bag