


Anti-Pedf/Serpinf1 Rabbit Monoclonal Antibody For Wb, Elisa - P36955

Item Number: CM0004328



Introduction

High-quality rabbit monoclonal antibody targeting human PEDF/SERPINF1 (P36955), validated for WB and ELISA. Detects PEDF in human, mouse, and rat samples with a molecular weight of 46kDa. Ideal for neurobiology and angiogenesis IVD studies.

[Learn More](#)

Parameter	Value
Product Name	PEDF/SERPINF1 Rabbit mAb
Remarks/Alias	O16; O112; PEDF; EPC-1; PIG35; PEDF/SERPINF1
Species	Human
Gene ID (Human)	5176
Gene ID	5176
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PEDF/SERPINF1 (P36955).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P36955
Protein Weight	46kDa
Shipping	Ice bag