


Anti-Pedf/Serpinf1 Rabbit Polyclonal Antibody For Wb, Elisa - P36955

Item Number: CM0014775



Introduction

Rabbit polyclonal antibody targeting human PEDF/SERPINF1, validated for WB and ELISA applications. Cross-reacts with mouse and rat. Suitable for angiogenesis and neurotrophic factor research.

[Learn More](#)

Parameter	Value
Product Name	PEDF/SERPINF1 Rabbit pAb
Remarks/Alias	OI6; OI12; PEDF; EPC-1; PIG35; PEDF/SERPINF1
Species	Human
Gene ID	5176
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 141-240 of human PEDF/SERPINF1 (NP_002606.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P36955
Protein Weight	46kDa
Shipping	Ice bag