


Anti-Pef1 Rabbit Monoclonal Antibody For Wb - Q9Ubv8

Item Number: CM0002862



Introduction

Recombinant Rabbit anti-PEF1 monoclonal antibody validated for WB and ELISA. Reacts with human, mouse, and rat PEF1. Ideal for studying ER-Golgi transport and calcium signaling.

[Learn More](#)

Parameter	Value
Product Name	PEF1 Rabbit mAb
Remarks/Alias	ABP32; PEF1A; PEF1
Species	Human
GeneID	553115
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human PEF1 (Q9UBV8)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBV8
Protein Weight	30kDa
Shipping	Ice bag