


Anti-Peg3 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Gzu2

Item Number: CM0016955



Introduction

PEG3 Rabbit pAb validates in WB, IHC-P, ELISA across human, mouse, rat. Recognizes PEG3 (181kDa), a zinc finger protein regulating apoptosis and tumor suppression. Ideal for neurobiology and cancer studies.

[Learn More](#)

Parameter	Value
Product Name	PEG3 Rabbit pAb
Remarks/Alias	PW1; ZNF904; ZSCAN24; ZKSCAN22; PEG3
Species	Human
Gene ID (Human)	5178
Gene ID	5178
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9GZU2
Protein Weight	181kDa
Shipping	Ice bag