

Anti-Peli3 Polyclonal Antibody For Wb, Elisa - Q8N2H9

Item Number: CM0027323

Introduction

Rabbit polyclonal antibody targeting human PELI3, suitable for Western blot (WB) and ELISA applications. PELI3 is an E3 ubiquitin ligase involved in innate immunity and TLR/IL-1 signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	PELI3 Rabbit pAb
Remarks/Alias	PELI3
Species	Human
GeneID	246330
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human PELI3 (NP_659502.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8N2H9
Protein Weight	51kDa
Shipping	Ice bag