


Anti-Pelp1 Polyclonal Antibody For Ihc-P, Elisa - Q8Izl8

Item Number: CM0015828



Introduction

PELP1 rabbit polyclonal antibody validated for ELISA and IHC-P, with cross-reactivity to human, mouse, and rat. Ideal for estrogen receptor signaling, cancer biology, and rRNA maturation studies.

[Learn More](#)

Parameter	Value
Product Name	PELP1 Rabbit pAb
Remarks/Alias	MNAR; P160; PELP1
Species	Human
Gene ID	27043
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 884-1130 of human PELP1 (Q8IZL8)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IZL8
Protein Weight	120kDa
Shipping	Ice bag