


## Anti-Pelp1 Rabbit Monoclonal Antibody For Wb And Elisa - Q8IZL8

Item Number: CM0009047



### Introduction

High-quality rabbit monoclonal antibody targeting human PELP1 (MNAR), validated for WB and ELISA. Cross-reacts with mouse and rat. PELP1 is a transcriptional coactivator involved in estrogen receptor signaling and cancer biology. Ideal for research use.

[Learn More](#)

| Parameter        | Value                                                                                                                    |
|------------------|--------------------------------------------------------------------------------------------------------------------------|
| Product Name     | PELP1 Rabbit mAb                                                                                                         |
| Remarks/Alias    | MNAR; P160; PELP1                                                                                                        |
| Species          | Human                                                                                                                    |
| GeneID (Human)   | 27043                                                                                                                    |
| GeneID           | 27043                                                                                                                    |
| Immunogen        | Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1000-1100 of human PELP1 (Q8IZL8). |
| Source           | Rabbit                                                                                                                   |
| Category         | Monoclonal Antibodies                                                                                                    |
| Application      | WB, ELISA                                                                                                                |
| Cross Reactivity | Human, Mouse, Rat                                                                                                        |
| SWISS            | Q8IZL8                                                                                                                   |
| Protein Weight   | 120kDa                                                                                                                   |
| Shipping         | Ice bag                                                                                                                  |