

Anti-Pelp1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - C9Jfv4

Item Number: CM0027100

Introduction

High-quality rabbit polyclonal antibody against human PELP1 (MNAR), validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Ideal for studying nuclear receptor signaling and transcriptional regulation.

[Learn More](#)

Parameter	Value
Product Name	PELP1 Rabbit pAb
Remarks/Alias	MNAR; P160; PELP1
Species	Human
Gene ID (Human)	27043
Gene ID	27043
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1021-1180 of human PELP1 (NP_055204.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	C9JFV4
Protein Weight	125kDa
Shipping	Ice bag