


Anti-Pen2/Psenen Rabbit Monoclonal Antibody For Wb, If/Icc, If-P, Elisa - Q9Nz42

Item Number: CM0009418



Introduction

Recombinant rabbit monoclonal antibody targeting human PEN2/PSENEN (UniProt Q9NZ42). Validated for WB, IF/ICC, IF-P, ELISA with cross-reactivity to mouse and rat. Ideal for gamma-secretase complex and Alzheimer's research.

[Learn More](#)

Parameter	Value
Product Name	PEN2/PSENEN Rabbit mAb
Remarks/Alias	PEN2; PEN-2; MDS033; ACNINV2; MSTP064; PEN2/PSENEN
Species	Human
Gene ID (Human)	55851
Gene ID	55851
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-101 of human PEN2/PSENEN (NP_758844.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NZ42
Protein Weight	12kDa
Shipping	Ice bag