


Anti-Pen2/Psenen Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nz42

Item Number: CM0014154



Introduction

Rabbit polyclonal antibody against PEN2/PSENEN, a gamma-secretase subunit. Validated for WB and ELISA; cross-reacts with human, mouse, rat. Ideal for Alzheimer's disease and Notch signaling research.

[Learn More](#)

Parameter	Value
Product Name	PEN2/PSENEN Rabbit pAb
Remarks/Alias	PEN2; PEN-2; MDS033; ACNINV2; MSTP064; PEN2/PSENEN
Species	Human
Gene ID (Human)	55851
Gene ID	55851
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-101 of human PEN2/PSENEN (NP_758844.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NZ42
Protein Weight	12 kDa
Shipping	Ice bag