

## Anti-Penk Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P01210

Item Number: CM0029017



### Introduction

Rabbit polyclonal antibody against human Proenkephalin-A (PENK), validated for WB, IF-P, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Ideal for opioid signaling and pain research.

[Learn More](#)

| Parameter        | Value                                                                                                                                    |
|------------------|------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | PENK Rabbit pAb                                                                                                                          |
| Remarks/Alias    | PE; PENK-A; PENK                                                                                                                         |
| Species          | Human                                                                                                                                    |
| GeneID (Human)   | 5179                                                                                                                                     |
| GeneID           | 5179                                                                                                                                     |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 25-267 of human PENK (NP_001129162.1). |
| Source           | Rabbit                                                                                                                                   |
| Category         | Polyclonal Antibodies                                                                                                                    |
| Application      | WB, IF-P, IHC-P, ELISA                                                                                                                   |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                        |
| SWISS            | P01210                                                                                                                                   |
| Protein Weight   | 31kDa                                                                                                                                    |
| Shipping         | Ice bag                                                                                                                                  |