


Anti-Penk Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P01210

Item Number: CM0022631



Introduction

Rabbit polyclonal antibody targeting human proenkephalin-A (PENK) for Western blot, immunofluorescence, and ELISA. Cross-reacts with rat. Validated for research on pain perception, opioid signaling, and neuropeptide function.

[Learn More](#)

Parameter	Value
Product Name	PENK Rabbit pAb
Remarks/Alias	PE; PENK-A; Proenkephalin-A
Species	Human
Gene ID (Human)	5179
Gene ID	5179
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 150-267 of human PENK (NP_001129162.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P01210
Protein Weight	31kDa
Shipping	Ice bag