

Anti-Pepsinogen I (Pg I) Rabbit Monoclonal Antibody For Elisa, Ihc-P, If-P - P0Djd7

Item Number: CM0007668



Introduction

Rabbit monoclonal antibody against human Pepsinogen I (PG I/PGA4) validated for ELISA, IHC-P, and IF-P. Cross-reacts with mouse and rat. Ideal for digestive enzyme and gastric cancer research.

[Learn More](#)

Parameter	Value
Product Name	Pepsinogen I (PG I) Rabbit mAb
Remarks/Alias	Pepsinogen-4; Pepsinogen I (PG I)
Species	Human
GeneID (Human)	643847
GeneID	643847
Immunogen	Recombinant protein of human Pepsinogen I (PG I)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P0DJD7
Protein Weight	42kDa
Shipping	Ice bag