

Anti-Pepsinogen Ii (Pg Ii) Rabbit Monoclonal Antibody For Ihc-P, Elisa - P20142

Item Number: CM0006595



Introduction

Rabbit monoclonal antibody against human Pepsinogen II (PG II), validated for IHC-P and ELISA. Cross-reacts with mouse. Recombinant immunogen. Suitable for digestive enzyme research and gastric cancer biomarker studies.

[Learn More](#)

Parameter	Value
Product Name	Pepsinogen II (PG II) Rabbit mAb
Remarks/Alias	PEPC; PGII; Pepsinogen II (PG II)
Species	Human
GeneID (Human)	5225
GeneID	5225
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-388 of human Pepsinogen II (PG II) (NP_002621.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P20142
Protein Weight	42kDa
Shipping	Ice bag