

Anti-Pepsinogen II (Pg II) Rabbit Monoclonal Antibody For IHC-P, Elisa - P20142

Item Number: CM0003084



Introduction

Rabbit monoclonal antibody against human Pepsinogen II (PG II). Suitable for IHC-P and ELISA, with cross-reactivity to mouse and rat. Ideal for gastrointestinal research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	Pepsinogen II (PG II) Rabbit mAb
Remarks/Alias	PEPC; PGII; Pepsinogen II (PG II)
Species	Human
GeneID (human)	5225
GeneID	5225
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-388 of human Pepsinogen II (PG II) (NP_002621.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20142
Protein Weight	42kDa
Shipping	Ice bag