

Anti-Per1 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - O15534

Item Number: CM0005556

Introduction

Anti-PER1 rabbit monoclonal antibody validated for WB, IF/ICC, ELISA. Targets human PER1 (O15534), a core circadian clock protein. Suitable for circadian rhythm and gene expression studies.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] PER1 Rabbit mAb
Remarks/Alias	PER; hPER; RIGUI
Species	Human
GeneID (Human)	5187
GeneID	5187
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 207-275 of human PER1 (NP_002607.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	O15534
Protein Weight	136 kDa/92 kDa
Shipping	Ice bag