

Anti-Per1 Rabbit Polyclonal Antibody - O15534

Item Number: CM0027056

Introduction

Rabbit polyclonal antibody against human PER1, validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. For circadian clock research and transcriptional regulation studies.

[Learn More](#)

Parameter	Value
Product Name	PER1 Rabbit pAb
Remarks/Alias	PER; hPER; RIGUI; PER1
Species	Human
Gene ID (Human)	5187
Gene ID	5187
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human PER1 (NP_002607.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15534
Protein Weight	136 kDa/92 kDa
Shipping	Ice bag