

Anti-Per3 Rabbit Monoclonal Antibody For Wb, Elisa - P56645

Item Number: CM0007548

Introduction

Rabbit monoclonal antibody targeting human PER3 circadian clock protein. Validated for Western blot and ELISA, cross-reacts with human, mouse, and rat. Ideal for sleep and circadian rhythm research.

[Learn More](#)

Parameter	Value
Product Name	PER3 Rabbit mAb
Remarks/Alias	GIG13; FASPS3; PER3
Species	Human
Gene ID (Human)	8863
Gene ID	8863
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 10-100 of human PER3 (P56645).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P56645
Protein Weight	132kDa
Shipping	Ice bag