

Anti-Perforin Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P14222

Item Number: CM0011017

Introduction

High-quality rabbit monoclonal antibody targeting human perforin (PRF1, P14222) for WB, IF/ICC, IHC-P, ELISA. Cross-reacts with mouse. Suitable for immunology and cancer research.

[Learn More](#)

Parameter	Specification
Product Name	Perforin Rabbit mAb
Remarks/Alias	P1, PFP, HPLH2
Species	Human
GenelD (human)	5551
GenelD	5551
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 21-555 of human Perforin (NP_005032.2))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P14222
Protein Weight	61 kDa
Shipping	Ice bag