

Anti-Perforin Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P14222

Item Number: CM0018649

Introduction

Rabbit polyclonal antibody against human Perforin (PRF1, P14222). Validated for WB, IHC-P, IF/ICC, ELISA with cross-reactivity to Human, Mouse, and Rat. A key tool for studying granzyme-mediated cytotoxicity and immune defense.

[Learn More](#)

Parameter	Value
Product Name	Perforin Rabbit pAb
Remarks/Alias	P1; PFP; HPLH2; Perforin
Species	Human
Gene ID (Human)	5551
Gene ID	5551
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 280-555 of human Perforin (NP_005032.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14222
Protein Weight	61 kDa
Shipping	Ice bag