

Anti-Perilipin-2 Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q99541

Item Number: CM0030218

Introduction

Rabbit polyclonal antibody against Perilipin-2 (ADRP/ADFP). Validated for WB, IF/ICC, IHC-P, ELISA. Reacts with human, mouse, rat. Lipid droplet structural component. Research use for lipid metabolism studies. UniProt Q99541.

[Learn More](#)

Parameter	Value
Product Name	Perilipin-2 Rabbit pAb
Remarks/Alias	ADFP; ADRP; Perilipin-2
Species	Human
GeneID (Human)	123
Immunogen	Recombinant protein A synthetic peptide corresponding to a sequence within amino acids 168-437 of human Perilipin-2 (NP_001113.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99541
Protein Weight	48 kDa
Shipping	Ice bag