

Anti-Perilipin-2 Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P43883

Item Number: CM0021535

Introduction

Rabbit polyclonal antibody targeting mouse Perilipin-2 (ADRP/Adipophilin) for WB, IF-P, IHC-P, and ELISA; cross-reacts with rat. Ideal for lipid droplet biology and metabolic disease research.

[Learn More](#)

Parameter	Value
Product Name	perilipin-2 Rabbit pAb
Remarks/Alias	ADPH, Adfp, Adrp, perilipin-2
Species	Mouse
GeneID (human)	123
GeneID	11520
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 325-425 of mouse perilipin-2 (NP_031434.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P43883
Protein Weight	46kDa
Shipping	Ice bag