

Anti-Perilipin A Monoclonal Antibody For Wb, Ihc-P, If-P, Elisa - O60240

Item Number: CM0009968

Introduction

Perilipin A Rabbit Monoclonal Antibody validated for WB, IHC-P, IF-P, and ELISA with cross-reactivity to mouse and rat. Ideal for adipocyte lipid metabolism research, targeting O60240.

[Learn More](#)

Parameter	Value
Product Name	Perilipin A Rabbit mAb
Remarks/Alias	PERI; PLIN; FPLD4; Perilipin A
Species	Human
GeneID (Human)	5346
GeneID	5346
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 400-500 of human Perilipin A (O60240)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O60240
Protein Weight	56 kDa
Shipping	Ice bag