


Anti-Periostin Monoclonal Antibody For Wb, Elisa - Q15063

Item Number: CM0011541



Introduction

Human Periostin rabbit monoclonal antibody validated for Western blotting and ELISA. Raised against a synthetic peptide within the N-terminus. Ideal for studying cell adhesion, tumor microenvironment, and extracellular matrix remodeling.

[Learn More](#)

Parameter	Value
Product Name	Periostin Rabbit mAb
Remarks/Alias	PN; OSF2; OSF-2; PDLPOSTN; Periostin
Species	Human
GeneID (Human)	10631
GeneID	10631
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Periostin (NP_001129406.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q15063
Protein Weight	80kDa/83kDa/84kDa/87kDa/90kDa/93kDa
Shipping	Ice bag