

Anti-Perk Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q9Nzj5

Item Number: CM0012542



Introduction

PERK Rabbit mAb targets human PERK (EIF2AK3), a key UPR sensor. Rabbit monoclonal, raised against amino acids 867-1116. Applicable for WB, IF/ICC, IHC-P, ELISA; cross-reacts with mouse and rat. For ER stress research.

[Learn More](#)

Parameter	Value
Product Name	PERK Rabbit mAb
Remarks/Alias	PEK; WRS; PERK
Species	Human
GeneID (Human)	9451
GeneID	9451
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 867-1116 of human PERK (NP_004827.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NZJ5
Protein Weight	125kDa
Shipping	Ice bag