

Anti-Perk Rabbit Monoclonal Antibody For Wb, Elisa - Q9Nzj5

Item Number: CM0008831



Introduction

PERK Rabbit Monoclonal Antibody (Q9NZJ5) validated for Western blot and ELISA, detecting human PERK with cross-reactivity to mouse and rat. Ideal for unfolded protein response and stress signaling research.

[Learn More](#)

Parameter	Value
Product Name	PERK Rabbit mAb
Remarks/Alias	PEK; WRS; PERK
Species	Human
GeneID (Human)	9451
GeneID	9451
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NZJ5
Protein Weight	125kDa
Shipping	Ice bag