


Anti-Peroxiredoxin 1/Pag Polyclonal Antibody For Wb And Elisa - Q06830

Item Number: CM0016387



Introduction

CamelBio's Peroxiredoxin 1/PAG Rabbit pAb (Q06830) is a polyclonal antibody validated for WB and ELISA in human and mouse samples. It detects the 22kDa antioxidant protein Peroxiredoxin-1, a key regulator of oxidative stress signaling.

[Learn More](#)

Parameter	Value
Product Name	Peroxiredoxin 1/PAG Rabbit pAb
Remarks/Alias	PAG; PAGA; PAGB; PRX1; PRXI; MSP23; NKEFA; TDPX2; NKEF-A; Peroxiredoxin 1/PAG
Species	Human
GeneID (Human)	5052
GeneID	5052
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-199 of human Peroxiredoxin 1/PAG (NP_001189360.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q06830
Protein Weight	22kDa
Shipping	Ice bag