


Anti-Peroxiredoxin 1/Pag Rabbit Mab For Wb, Elisa - Q06830

Item Number: CM0008689



Introduction

Rabbit monoclonal antibody against Peroxiredoxin 1/PAG (PRDX1, NKEF-A). Detects human and mouse PRDX1 for Western blot and ELISA applications. Useful for oxidative stress and antioxidant research.

[Learn More](#)

Parameter	Value
Product Name	Peroxiredoxin 1/PAG Rabbit mAb
Remarks/Alias	PAG; PAGA; PAGB; PRX1; PRXI; MSP23; NKEFA; TDPX2; NKEF-A; Peroxiredoxin 1/PAG
Species	Human
Gene ID (Human)	5052
Gene ID	5052
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Peroxiredoxin 1/PAG (Q06830)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q06830
Protein Weight	22kDa
Shipping	Ice bag