


Anti-Peroxiredoxin 1/Pag Rabbit Polyclonal Antibody For If/Icc, Elisa - Q06830

Item Number: CM0020474



Introduction

Rabbit polyclonal antibody targeting Peroxiredoxin 1/PAG for IF/ICC and ELISA applications. Cross-reacts with human, mouse, and rat. Detects a 22kDa protein. Suitable for oxidative stress and signaling research.

[Learn More](#)

Parameter	Value
Product Name	Peroxiredoxin 1/PAG Rabbit pAb
Remarks/Alias	PAG; PAGA; PAGB; PRX1; PRXI; MSP23; NKEFA; TDPX2; NKEF-A; Peroxiredoxin 1/PAG
Species	Human
GeneID (Human)	5052
GeneID	5052
Immunogen	Recombinant protein A synthetic peptide corresponding to a sequence within amino acids 99-199 of human Peroxiredoxin 1/PAG (NP_001189360.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q06830
Protein Weight	22kDa
Shipping	Ice bag