


Anti-Peroxiredoxin 2 (Prdx2) Monoclonal Antibody For Wb And Elisa - P32119

Item Number: CM0007917



Introduction

High-quality rabbit monoclonal antibody against human Peroxiredoxin 2 (PRDX2), validated for Western blot and ELISA. Reacts with human, mouse, and rat. Ideal for oxidative stress and redox research.

[Learn More](#)

Parameter	Value
Product Name	Peroxiredoxin 2 (PRDX2) Rabbit mAb
Remarks/Alias	PRP; TSA; PRX2; PTX1; TPX1; NKEFB; PRXII; TDPX1; NKEF-B; HEL-S-2a; Peroxiredoxin 2 (PRDX2)
Species	Human
GeneID (Human)	7001
GeneID	7001
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P32119
Protein Weight	22kDa
Shipping	Ice bag