


Anti-Peroxiredoxin 2 (Prdx2) Polyclonal Antibody For Wb, If/Icc, Elisa - P32119

Item Number: CM0022034



Introduction

Rabbit polyclonal antibody targeting Peroxiredoxin 2 (PRDX2), validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat; immunogen recombinant protein. Ideal for oxidative stress and antioxidant research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Peroxiredoxin 2 (PRDX2) Rabbit pAb
Remarks/Alias	PRP; TSA; PRX2; PTX1; TPX1; NKEFB; PRXII; TDPX1; NKEF-B; HEL-S-2a; 2)
Species	Human
Gene ID	7001
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-142 of human Peroxiredoxin 2 (PRDX2) (NP_859428.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P32119
Protein Weight	22kDa
Shipping	Ice bag