


Anti-Peroxiredoxin 2 (Prdx2) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P32119

Item Number: CM0024030



Introduction

Anti-PRDX2 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA in human, mouse, rat. Detects ~22 kDa peroxiredoxin-2, a cytoplasmic antioxidant enzyme involved in oxidative stress response.

[Learn More](#)

Parameter	Value
Product Name	Peroxiredoxin 2 (PRDX2) Rabbit pAb
Remarks/Alias	PRP; TSA; PRX2; PTX1; TPX1; NKEFB; PRXII; TDPX1; NKEF-B; HEL-S-2a; Peroxiredoxin 2 (PRDX2)
Species	Human
Gene ID	7001 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-142 of human Peroxiredoxin 2 (PRDX2) (NP_859428.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P32119
Protein Weight	22kDa
Shipping	Ice bag