


Anti-Peroxiredoxin 3 (Prdx3) Polyclonal Antibody For Wb, If, Elisa - P30048

Item Number: CM0014729



Introduction

Rabbit polyclonal antibody targeting human PRDX3 (Peroxiredoxin 3), validated for WB, IF/ICC, and ELISA. Detects endogenous PRDX3 across human, mouse, and rat samples. Ideal for oxidative stress and mitochondrial research.

[Learn More](#)

Parameter	Value
Product Name	Peroxiredoxin 3 (PRDX3) Rabbit pAb
Remarks/Alias	AOP1; MER5; AOP-1; PPPCD; SP-22; HBC189; SCAR32; PRO1748; prx-III; Peroxiredoxin 3 (PRDX3)
Species	Human
GeneID (Human)	10935
GeneID	10935
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 63-256 of human Peroxiredoxin 3 (Peroxiredoxin 3 (PRDX3)) (NP_006784.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30048
Protein Weight	28kDa
Shipping	Ice bag